



- \* In the outdoors:
- \* During sporting activities:
- \* When toilets are unclear:
- \* To avoid huge queues:
- \* By the mobility impaired:
- \* In every day life:
- \* In clinical situations:

### When is a Whiz® used?

for when nature calls



FOR FURTHER INFORMATION  
SEE OVERLEAF OR VISIT  
[www.whizfreedom.com](http://www.whizfreedom.com)



Authorised Representative:  
Freedomatlast Ltd, Park Manor,  
Upper Malrow St. Limerick, Ireland  
[www.whizaway.com](mailto:info@whizaway.com) - [info@whizaway.com](mailto:info@whizaway.com)  
CE marked and FDA approved  
Available on NHS prescription  
Whiz is a registered trademark of Jbol Ltd.

Manufactured & Distributed by: JBOL Ltd  
1 Folly Bridge  
Oxford OX1 4LB, UK  
Tel: +44 (0) 1865 240572  
Fax: +44 (0) 1865 246565  
[www.jbol.co.uk](http://www.jbol.co.uk)  
[info@jbol.co.uk](mailto:info@jbol.co.uk)

+ For all medical uses, including incontinence, see website or contact us.

- \* Camping, hiking, hill walking, mountaineering
- \* Sailing, skiing, canoeing, scuba diving, running, biking
- \* Events, public loos, clubbing, gas stations, travelling
- \* Festivals, stadiums, public events, outdoor shows, public buildings
- \* When toilets are unavailable or difficult to use
- \* Traffic jams, on holiday, when you simply get caught short
- \* During pregnancy, arthritis, kidney evaluation.†



Giving all women the freedom to wee standing, sitting or crouching, indoors or outdoors, at their convenience without the need to remove clothing.

[www.whizfreedom.com](http://www.whizfreedom.com)

### Every woman should have one™

Used by women all over the world for convenience, travel, sport or medical needs. The perfect accessory for backpacking, festivals and a wide range of outdoor and sporting activities.

- Hygienic, antibacterial, latex-free, medical grade material
- Soft and comfortable, easy to use urine director
- Liquid repellent, flick to dry
- Lightweight, flexible, easy to store and carry
- Re-usable, machine washable and recyclable





Patents Granted: GB2362577  
GB2396819, GB23855276  
UK Reg. Design: 2100623, 2103229  
International Patent Application:  
01/74275, PCT/GB2003/000472  
Further patents pending © JBOL Ltd. Nov 2009

## Instructions for Use

Thank you for buying the Whiz<sup>®</sup> freedom.



1. Hold the Whiz<sup>®</sup> snugly against the body as shown



2-3. Best positions when using the Whiz<sup>®</sup>



4. Using the Whiz<sup>®</sup> without undressing

**Important:** Always position the Whiz<sup>®</sup> with the broad area of the shield pointing upwards. Held correctly the Whiz<sup>®</sup> will not leak and stays completely dry on the exterior. Keep the nozzle pointed downwards and away from the body. Do not share a Whiz<sup>®</sup> with others.

**Uses:** For the hygienic voiding and disposal of urine as a urine director; for women with continence issues; for renal or transplant patients; for men with penile amputation or retraction; for transgenders; for women with mobility issues or during pregnancy; and for when wcs are dirty or unavailable.

For full list of uses/users visit: [www.whizfreedom.com](http://www.whizfreedom.com)

**(4) Medical needs.** If you are bed bound attach a connector tube and either lie on your back and elevate your pelvis with a cushion, or roll onto your side and tuck your upper leg closer to your torso than your bottom leg. Hold the Whiz<sup>®</sup> in place and make sure the Whiz<sup>®</sup> and tube slant downwards. If you are in a chair or a wheelchair edge forward to the front of the seat with the Whiz<sup>®</sup> (and tube if applicable) slanting downwards.

For free downloadable tips, full medical IFU, brochures, other uses and testimonials from users, visit [www.whizfreedom.com](http://www.whizfreedom.com)



Hold your Whiz<sup>®</sup> this way up



**Tips on using the Whiz<sup>®</sup>:-**

**(1) Find your most comfortable position:** The Whiz<sup>®</sup> is easy to use, but still try it for the first time in the privacy of your home (in a shower for example) to find the best and most comfortable position for you – standing, sitting or crouching.

**(2) No need to undress:** the Whiz<sup>®</sup> can be used with a skirt, or trousers with a zip, without the need to undress, or lower clothing to the ankle. Simply pull your underwear to one side, or slightly lower it, and then wee normally through the Whiz<sup>®</sup> aiming away from the body.

**(3) Attaching a connector tube:** (an accessory from [whizfreedom.com](http://whizfreedom.com)): Useful if you are in a confined space, mobility impaired, or need to use a drainage bag. Push the end of the Whiz<sup>®</sup> nozzle inside the connector tube. It will be a tight fit. After usage shake the connector tube before storage. Remember to always position the end of the connector tube and drainage bag/bottle lower than the Whiz<sup>®</sup>.

**Caring for your Whiz<sup>®</sup>:-**

**(i) Replacement:** The Whiz<sup>®</sup> is liquid repellent, washable and re-usable (by a single person). Given the intimate nature of the Whiz<sup>®</sup> we recommend replacing it every

year, or after extensive use. NOTE: urine is sterile when it exits the body and poses no health hazard to the user.

**(ii) Cleaning:** After use, one flick and the Whiz<sup>®</sup> will be dry. Once home, the Whiz<sup>®</sup> should be rinsed in water. It can also be washed with regular washing up liquid, or in a washing machine (up to 40°C).

**(iii) Storage:** The Whiz<sup>®</sup> is light-weight and flexible. It can easily be rolled, or folded. The material will spring back to its original shape as long as it is not kept distorted for long periods. Always store your Whiz<sup>®</sup> unrolled. Bags and other accessories are available from [www.whizfreedom.com](http://www.whizfreedom.com)

**Updates:** Enjoy your Whiz<sup>®</sup> and feel free to give us your feedback. Product and accessories updates are available at [www.whizfreedom.com](http://www.whizfreedom.com)

**Technical:** The Whiz<sup>®</sup> is made from a medical-grade thermoplastic elastomer, which includes an anti-bacterial agent. It is hydrophobic (liquid repellent) and latex free. It is CE marked, FDA approved and is available on NHS prescription in the UK.

The Whiz<sup>®</sup> is patent and design right granted and protected. The manufacturer JBOL Ltd, have taken every care and attention to ensure that the Whiz<sup>®</sup> freedom works as effectively as possible. However, under no circumstances will JBOL Ltd be held liable for any damage, injury or inability to use the Whiz<sup>®</sup> freedom.

